

REASONS FOR USING THE EZ-CLEAN

Colostomy / Ileostomy POUCH SYSTEM

1. The EZ-Clean System closely duplicates normal body function disposing of fecal matter in the toilet to be flushed away.
2. The EZ-Clean System helps the person adjust to having an ostomy and the mental strain of not having control of bowel function and waste disposal.
3. It is fast, hygienic and the user does not have to come in contact with fecal matter.
4. The pouch was designed to duplicate an individual's normal bowel movement function and waste disposal.
5. The EZ-Clean pouch was designed by an ostomate to normalize his life after undergoing colorectal surgery that resulted in a permanent Colostomy.
6. The EZ-Clean pouch was engineered with an internal dispersion manifold that sprays water under pressure from the top down to thoroughly cleanse the pouch
7. Waste is flushed out of the pouch drain into the toilet.
8. The entire process takes under 3 minutes and the user wipes only clean water from the drain before closing the pouch.
9. The pouch water inlet also serves to allow the release of accumulated gas in the pouch when necessary.
10. The EZ-Clean pouch is one piece, very flexible and only 3/16 inch thick.
11. There are no bulky couplings to snap together creating bulges under clothing and obstructions when you try to bend over.
12. The ostomate does not have to squeeze out body waste and spend inordinate amounts of time in the bathroom trying to clean out the pouch.
13. The pouch can be worn for all daily activities including; showers, baths, work, athletic events for 3 to 5 days or longer before replacement.
14. No special skin preps are required for the EZ-Clean System, only clean dry abdominal skin to allow a secure attachment.
15. The user sits on the toilet, opens the drain, attaches the system and flushes the pouch.
16. All fecal matter is deposited in the toilet and flushed away as would occur with a normal person after a bowel movement.
17. No bagging pouches full of fecal matter for disposal in the trash (They are not flushable in the toilet).

18. No storing the weeks bagged pouches full of human waste for trash collection.
19. No rinsing fecal matter in the sink from two piece pouch systems or hanging them in the bathroom to dry for exchange when the one being worn is full.
20. No household odor from pouches full of fecal matter stored in plastic bags
21. No mental apprehension about spending hours each day in the bathroom trying to clean waste from the pouch.
22. Peace of mind knowing the cleaning process will not be an unknown ordeal without any idea how long it will take.
23. Having a choice of three (3) different water feed systems to choose from to deal with any situation:
 - A. Fixed Feed Water System (Hooks up to the water closet inlet of the toilet and hangs ready for use when needed).
 - B. Portable Feed Water System (Attaches to the faucet of a sink in place of the screen used to collect foreign material in the water) used for travel and is adaptable to hotel / motel stays.
 - C. Water Bottle Adapter (Plugs into most 6 – 8 ounce plastic water bottles to permit cleaning at any toilet facility). An aspiration valve prevents siphoning back pouch contents.
24. A well thought out system designed to bring normality back to your life after undergoing life changing surgery.
25. Revolutionary design that was awarded the 21 Century Health Award in 2012 for its unique carefully thought-out engineering.
26. Why would anyone choose an obsolete system to care for their Colostomy / Ileostomy when they can have a modern product that normalizes their life?
27. Why take a bumpy road when you can use a smooth comfortable road for travel.
28. There are company subsidized “Trial Kits” available which permit ostomates to experience the miracle of this life changing product.
29. More new products are coming from this manufacturer to make life easier for ostomates.
30. The product is protected by 5 US and 4 Foreign patents and is the only system of its kind; available only from Schena Ostomy Technologies or one of the appointed medical distributors who carry their products.
31. A new patented dispersion manifold is now available which will allow ostomates to modify their drainable pouch for use with the EZ-Clean System.

What are you waiting for? Make your life easier and normal!

